

# TUBE SEALER

## Model SE700



### Features

- Battery-powered, compact and lightweight.
- Seals at handheld sealing unit.
- Sealing unit cable is 1.8 meters (6 feet) long.
- Seals wide range of blood bag tubing - good for tubing diameter 2-6 mm (1/4").
- Fast sealing in one second. No warm-up required.
- Makes seals with a notch at the center for easy segmenting of tubing.
- Full charging in 6 hours. Battery level meter indicates charge level and charging status.
- Battery pack is rechargeable and replaceable. Battery charger is built-in in the unit.
- Head protector is detachable for sealing head cleaning.
- Canvas carrying bag is included. The bag carries all parts.
- Lightweight, battery operation and carrying bag - an ideal device for mobile use.

## Specification

- (1) Sealing time : 1 - 2 second
- (2) Tubing diameter : 2 - 6 mm (1/4")
- (3) RF cable length : 1.8 m (6')
- (4) Indication lamp : BATTERY LEVEL, CHARGING, SEAL, SEALING
- (5) Operation frequency : 40.68 MHz
- (6) Power source : battery pack (NiMH, 12VDC/4.5AH)
- (7) Battery charger unit : 100-120/220-240 VAC, 18 VDC/1 Amp out, 25 W
- (8) Battery charging : full charging in 6 hours
- (9) Dimension - main unit : W 148 x H 34 x D 252 mm (W 5.8" x H 1.4" x D 9.9")  
- hand unit : W 28 x H 44 x D 204 mm (W 1.1" x H 1.8" x D 8.1")
- (10) Weight - main unit : 1.7 Kg (3.8 lb)  
- hand unit : 0.2 Kg (0.5 lb)
- (11) Temperature - operation : 0 - 40 °C (32 - 104 °F)  
storage : -20 - 70 °C (-4 - 158 °F)

## Accessories



RF cable



Charging adapter



Carrying bag



Battery charge level is displayed on the panel.



The carrying bag houses all parts.